
# **Dental Stools**



#### D3 Doctor stool

Chrome plated pneumatic cylinder
Height adjustment,
Backrest forward backward adjustment,
Seat surface adjustment,
Contoured seat & back,
Fiber leather upholstery
Cast aluminum star,
Quiet twin wheel casters,
1 year warranty




Part number: A080103

Part number: A080224



#### **N4** Assistant stool

Round seat,
Fiber leather upholstery,
Chrome plated pneumatic cylinder
Adjustable chrome foot ring
Cast aluminum star
Ratcheting body support,
1 year warranty





Carrier Technical support: 800 488 9708 Add: 1590 S Milliken Ave., unit A, Ontario, CA 91761



**ADS Dental System** 





# **AJ16 Chair**

Part number: A091602

Luxury hydraulic chair with 3 programmable positions.

Reliable hydraulic system makes powerful and stable movement.

5 year warranty.



# AJ16 Chair

- Luxury Modern Superior
- With 3 programmable positions and reliable system,
- AJ16 moves stably, providing exceptional comfort for patients and helping them release.
- Swing up/down aluminum armrest firmly holds the patient and facilitates the patient to enter and exit.



 Wide backrest, cast aluminum backrest stand fully support the patient's back with ergonomic design.
 Butterfly shaped back offers wide space for doctors and enables you get as close to the patient as possible.



 Japanese hydraulic pump provides sufficient power to support the chair's movement.



Articulating headrest
 offers perfect support for the patient.
 Adjust the angle and get the best view
 of the oral cavity.



 60 degree seat rotations satisfy all kinds of needs.



 Cast aluminum base enables the chair moves steadily.



 Foot control makes operation more convenient for doctors.



## Specifications:

Environmental requirements — Ambient temperature 5 °C~40 °C;

Humidity≤80%;

Atomospherical pressure 8.7 psi~ 15.4psi;


Electrical requirements — AC110V ± 22V;


50HZ, 60HZ;


Hydraulic pump run time —— Intermittent operation 2min On / 18min Off;

Weight — 366lbs;

Max patient load — 550lbs;









Part number: A9161002



### Standard Features:

AJ16 Luxury hydraulic chair, with 3 programmable positions,
Articulating headrest, Swing up/down aluminium armrest,
Wide backrest, upgraded fiber leather, Cast iron chair frame,
Cast aluminium backrest stand, Cast aluminium base,
60 degrees seat rotation, chair safety switch,
One touchpad, Foot control

Classic dentist's control head with 3 auto H.P. control, Individual water coolant controls,

Master handpiece air coolant adjustment,

Asepsis handpiece tubing,

Master on/off switch, 1 doctor syringe, Wet/Dry foot control, Stainless steel tray with non-slip pad, Flex arm air brake switch,

Classic 100 Post mounted utility (PMU) with steel side box, City water/bottled system, Water outlet and air QD, Solids collector,

2 HVE, 1SE, 1 syringe on telescoping assistant's control,
Junction box with air/water filters, regulators, master shut-off valves,
5 ft. umbilical

Part number: A9161012



### **Standard Features:**

AJ16 Luxury hydraulic chair, with 3 programmable positions, Articulating headrest, Swing up/down aluminium armrest, Wide backrest, upgraded fiber leather, Cast iron chair frame, Cast aluminium backrest stand, Cast aluminium base, 60 degrees seat rotation, chair safety switch, One touchpad, Foot control

Classic dentist's control head with 3 auto H.P. control, Individual water coolant controls,

Master handpiece air coolant adjustment,

Asepsis handpiece tubing,

Master on/off switch, 1 doctor syringe,

Remote toggles for cup fill and bowl rinse, Wet/Dry foot control,

Stainless steel tray with non-slip pad, Flex arm air brake switch,

Porcelain cuspidor bowl with strainer, Cup fill & Bowl rinse toggles.

Classic 100 Post mounted utility (PMU) with steel side box,

City water/ bottled system, Water outlet and air QD,

Solids collector,

2 HVE, 1SE, 1 syringe on telescoping assistant's control, Junction box with air/water filters, regulators, master shut-off valves, 5 ft. umbilical

Part number: A9162002



## Standard Features:

AJ16 Luxury hydraulic chair, with 3 programmable positions, Articulating headrest, Swing up/down aluminium armrest, Wide backrest, upgraded fiber leather, Cast iron chair frame, Cast aluminium backrest stand, Cast aluminium base, 60 degrees seat rotation, chair safety switch,

One touchpad, Foot control

Swing Classic dentist's control head with 3 auto H.P. control, Individual water coolant controls,

Master handpiece air coolant adjustment,

Asepsis handpiece tubing,

Master on/off switch, 1 doctor syringe,

Wet/Dry foot control, Stainless steel tray with non-slip pad,

City/bottled water system, 5 ft umbilical,

Junction box with air/water filters, regulators, master shut-off valves,

Chair mounted swing rear assistant vacumm package,

Water outlet and air QD on assistant tray,

Solids collector,

2 HVE, 1SE, 1 syringe

Part number: A9162012



#### **Standard Features:**

AJ16 Luxury hydraulic chair, with 3 programmable positions, Articulating headrest, Swing up/down aluminium armrest, Wide backrest, upgraded fiber leather, Cast iron chair frame, Cast aluminium backrest stand, Cast aluminium base, 60 degrees seat rotation, chair safety switch,

One touchpad, Foot control

Swing Classic dentist's control head with 3 auto H.P. control, Individual water coolant controls,

Master handpiece air coolant adjustment,

Asepsis handpiece tubing,

Master on/off switch, 1 doctor syringe,

Wet/Dry foot control, Stainless steel tray with non-slip pad,

City/bottled water system, 5 ft umbilical,

Junction box with air/water filters, regulators, master shut-off valves,

Swing mount Cupspidor and assistant vacumm package,

Water outlet and air QD on assistant tray,

Solids collector,

2 HVE, 1SE, 1 syringe

Part number: A9163002



## Standard Features:

AJ16 Luxury hydraulic chair, with 3 programmable positions, Articulating headrest, Swing up/down aluminium armrest, Wide backrest, upgraded fiber leather, Cast iron chair frame, Cast aluminium backrest stand, Cast aluminium base, 60 degrees seat rotation, chair safety switch,

One touchpad, Foot control

Beyond luxury dentist's control head, Standard touchpad, Cast aluminium delivery head base, double pivot delivery head, with 3 auto H.P. control,

Individual water coolant controls,

Master handpiece air coolant adjustment,

Asepsis handpiece tubing,

Master on/off switch, 1 doctor syringe, Wet/Dry foot control, Stainless steel tray with non-slip pad, Flex arm air brake switch Beyond 300 Post mounted utility (PMU) with silver side box, City water/bottled system, Water outlet and air QD, Swing assistant tray with touchpad,

Solids collector,

2 HVE, 1SE, 1 syringe, Alumium base Junction box with air/water filters, regulators, master shut-off valves, 6' internal umbilical

Part number: A9163012



## **Standard Features:**

AJ16 Luxury hydraulic chair, with 3 programmable positions, Articulating headrest, Swing up/down aluminium armrest, Wide backrest, upgraded fiber leather, Cast iron chair frame, Cast aluminium backrest stand, Cast aluminium base, 60 degrees seat rotation, chair safety switch, One touchpad, Foot control

Beyond luxury dentist's control head, Standard touchpad, Cast aluminium delivery head base, double pivot delivery head, with 3 auto H.P. control,

Individual water coolant controls,

Master handpiece air coolant adjustment,

Asepsis handpiece tubing,

Master on/off switch, 1 doctor syringe, Wet/Dry foot control, Stainless steel tray with non-slip pad, Flex arm air brake switch, Porcelain cuspidor bowl with strainer, Cup fill & Bowl rinse toggles. Beyond 301 Post mounted utility (PMU) with steel side box, City water/bottled system, Water outlet and air QD, Solids collector,

2 HVE, 1SE, 1 syringe on telescoping assistant's control, Junction box with air/water filters, regulators, master shut-off valves, 5 ft. umbilical

Part number: A9164002



## Standard Features:

AJ16 Luxury hydraulic chair, with 3 programmable positions,
Articulating headrest, Swing up/down aluminium armrest,
Wide backrest, upgraded fiber leather, Cast iron chair frame,
Cast aluminium backrest stand, Cast aluminium base,
60 degrees seat rotation, chair safety switch,
One touchpad attached on armrest, Foot control
Swing Beyond dentist's control head, Standard touchpad,
Cast aluminium delivery head base, double pivot delivery head,
3 auto H.P. control, Individual water coolant controls,
Master handpiece air coolant adjustment, Asepsis handpiece tubing,

Master mandpiece air coolant adjustment, Asepsis handpiece tubing, Master on/off switch, 1 doctor syringe,

Wet/Dry foot control, Stainless steel tray with non-slip pad,

City/bottled water system, 5 ft umbilical,

Junction box with air/water filters, regulators, master shut-off valves,

Chair mounted swing rear assistant vacumm package,
Assistant touchpad, Water outlet and air QD on assistant tray,
Solids collector, 2 HVE, 1SE, 1 syringe
Post mount LED light with sensor control, 3 level intensity
adjustment

Part number: A9164012



## **Standard Features:**

AJ16 Luxury hydraulic chair, with 3 programmable positions, Articulating headrest, Swing up/down aluminium armrest, Wide backrest, upgraded fiber leather, Cast iron chair frame, Cast aluminium backrest stand, Cast aluminium base, 60 degrees seat rotation, chair safety switch,

One touchpad, Foot control

Swing Beyond dentist's control head, Standard touchpad, Cast aluminium delivery head base, double pivot delivery head, 3 auto H.P. control,

Individual water coolant controls,

Master handpiece air coolant adjustment,

Asepsis handpiece tubing,

Master on/off switch, 1 doctor syringe, Wet/Dry foot control,

Stainless steel tray with non-slip pad,

City/bottled water system, 5 ft umbilical,

Junction box with air/water filters, regulators, master shut-off valves,

Chair mounted swing rear cuspidor and assistant vacumm

Assistant touchpad, Water outlet and air QD on assistant tray, Solids collector,

2 HVE, 1SE, 1 syringe