



Amber II LED Light with arm and bushing
 Part number
A0600670



Post mount Amber II LED Dental Light 30"/40"
 Part number
 30" Post mount **A0602670**
 40" Ortho Post mount **A0602672**



**Cabinet mount Amber II LED Dental Light
 (cabinet not included)**
 Part number
A0601674



Wall mount Amber II LED Dental Light
 Part number
A0606670



Ceiling mount Amber II LED Dental Light
 Part number
A0601671 8' Ceiling mount
A0601672 9' Ceiling mount
A0601673 10' Ceiling mount

Amber II Dental LED Light



2023 version1



Amber II

Dental LED Light

Standard Features

Power input: AC24V
Illuminance: 8000-40000 lux

Color temperature
White light: 5200k
Yellow Light: 3500k

Composite Light: 4400k-5200k
Power consumption: <10W

Amber II LED Light head

Part number

A121715



3 axis head rotation



Composite light



Yellow light



White light

Easy switch of 3 operating modes through the upper sensor White light, yellow light and composite light



on/off control



Intensity control
Adjusting color temperature through touching the upper sensor